

Customer Bulletin

Patient Monitoring & Connectivity



Windows Service Packs and Security Patches

These are the verified and validated Windows service packs and security patches for ICS G2 (v4.01.02)

ICS Server, All Services, and Clinical Access Client

Database Server and Services

OS_Version_SP Windows Server 2003 R2 Standard Edition SP2
System_Type Server, Stand-alone, Terminal Server
Country United States

IE_Version 7.0.5730.13
SQL Version SP SQL 2005 SP3 (9.00.4035.00)

Symantec_Antivirus (none)
(none)
10.1.5000.5
(none)

Hotfixes KB924667-v2; KB925398_WMP64; KB925876; KB925902;
KB926122; KB927891; KB929123; KB930178; KB931784;
KB932168; KB933729; KB933854; KB935839; KB935840;
KB936021; KB936357; KB936782; KB938127; KB938127-IE7;
KB938464; KB941202; KB941568; KB941569; KB941644;
KB942763; KB942830; KB942831; KB942840; KB943055;
KB943460; KB943485; KB943729; KB944338-v2; KB944533;
KB944653; KB945553; KB946026; KB948496; KB950760;
KB950762; KB950974; KB951066; KB951698; KB951748;
KB952069; KB952954; KB954211; KB954550-v5; KB954600;
KB955069; KB955839; KB956802; KB956803; KB956841;
KB957097; KB958215; KB958644; KB958687; KB960714;
KB960715; KB961118; KB961260-IE7; KB967715; Q147222;

Clinical Access Client, Internet Explorer 7

OS_Version_SP Windows XP Professional SP3
System_Type Workstation
Country United States

IE_Version 7.0.5730.13
SQL Version SP (none)
Symantec_Antivirus (none)
(none)
10.1.5000.5

Customer Bulletin

Patient Monitoring & Connectivity



Hotfixes (none)
KB891122; KB909520; KB911564; KB925398_WMP64; KB929399;
KB936782_WMP11; KB936782_WMP9; KB938127-v2-IE7;
KB938464; KB939683; KB941569; KB943729; KB946648;
KB950759; KB950760; KB950762; KB950974; KB951066;
KB951376; KB951376-v2; KB951698; KB951748; KB951978;
KB952069_WM9; KB952287; KB952954; KB954154_WM11;
KB954211; KB954459; KB954550-v5; KB954600; KB955069;
KB955839; KB956391; KB956802; KB956803; KB956841;
KB957095; KB957097; KB958215; KB958215-IE7; KB958644;
KB958687; KB959772_WM11; KB960714; KB960714-IE7;
KB960715; KB961118; KB961260-IE7; KB967715; Q147222;
ServicePackUninstall;

OS_Version_SP © Windows Vista™ Enterprise SP1
System_Type Workstation, Terminal Server
Country United States

IE_Version 7.0.6001.18000
SQL Version SP (none)
Symantec_Antivirus (none)
(none)
(none)
(none)

Hotfixes KB937286, 943729, KB925902, KB929399, KB929451, KB929735,
KB929761, KB930163, KB930178, KB930857, KB931099,
KB931573, KB932394, KB932471, KB933579, KB933729,
KB936021, KB936357, KB936782, KB936825, KB937077,
KB938127, KB938952, KB939159, KB941202, KB941568,
KB941569, KB941600, KB943055, KB943078, KB945553,
KB946026, KB946456, KB905866, KB929123, KB929916,
KB931213, KB933928, KB935509, KB935807, KB936824,
KB937287, KB938123, KB938194, KB938371, KB938464,
KB938979, KB941649, KB941651, KB942624, KB942763,
KB943302, KB943411, KB943899, KB944533, KB946041,
KB948609, KB948610, KB950760, KB950762, KB950974,
KB951066, KB951376, KB951698, KB951978, KB952069,
KB952287, KB952709, KB953155, KB953733, KB954154,
KB954211, KB954459, KB955020, KB955069, KB955302,
KB955839, KB956802, KB956841, KB957097, KB957200,
KB957321, KB957388, KB958481, KB958483, KB958623,
KB958624, KB958644, KB958687, KB959108, KB959130,
KB959772, KB960715, KB961260, KB936330, 940157

Customer Bulletin

Patient Monitoring & Connectivity



Clinical Access Client, Internet Explorer 6

OS_Version_SP Windows XP Professional SP3
System_Type Workstation
Country United States

IE_Version 6.0.2900.5512;SP3;

SQL_Version_SP (none)

Symantec_Antivirus (none)

Hotfixes KB873339; KB885835; KB885836; KB886185; KB887472; KB888302;
KB890046; KB890859; KB891122; KB891781; KB893756; KB893803v2;
KB894391; KB896344; KB896358; KB896423; KB896428; KB898461;
KB899587; KB899591; KB900485; KB900725; KB901017; KB901214;
KB902400; KB904942; KB905414; KB905749; KB908519; KB908531;
KB909520; KB910437; KB911280; KB911562; KB911564; KB911927;
KB913580; KB914388; KB914389; KB916595; KB916846; KB917344;
KB918118; KB918439; KB919007; KB920213; KB920342; KB920670;
KB920683; KB920685; KB920872; KB922582; KB922819; KB923191;
KB923414; KB923980; KB924270; KB924496; KB924667;
KB925398_WMP64; KB925720; KB925876; KB925902; KB926239;
KB926255; KB926436; KB927779; KB927802; KB927891; KB928255;
KB928843; KB929123; KB929399; KB930178; KB930916; KB931261;
KB931784; KB932168; KB933729; KB935839; KB935840; KB936021;
KB936357; KB936782_WMP11; KB936782_WMP9; KB937894;
KB938127; KB938828; KB938829; KB939683; KB941202; KB941568;
KB941569; KB941644; KB941693; KB942763; KB942840; KB943055;
KB943460; KB943485; KB944338; KB944533; KB944653; KB945553;
KB946026; KB947864; KB948590; KB948881; KB950749; Q147222